

PARA PULSER®

TAKE BACK YOUR POWER



Product description

Ver 5.0

April 2026



PARAPULSER®

– the magnetic-(r)evolution

The **Parapulser®** is since 2012 the most compact, advanced and powerful magnetic pulser or **PEMF** device in its class (**P**ulsed **E**lectro **M**agnetic **F**ield). It is suitable for an almost unlimited range of applications in which strong magnetic pulses are used to activate the body's self-healing abilities. The **Parapulser®** is also particularly well suited for lymph node treatment according to Dr Beck's protocol.

The effect of the magnetic pulses is based on the physical law of '*magnetic induction*'. This causes all ions (i.e. charged atoms and molecules) within the body's cells to be set in motion. These so-called microcurrents increase the cell potential, thereby gently activating the body's self-healing powers.

To enhance the overall effect, the magnetic pulses generated by the **Parapulser®** are charged with informations from the quantum field using a golden spiral.



Parapulser®



Quantum field spiral

Technical highlights:

- **Pulse frequency 1 – 20 Hz** : With 12 intensity-steps in the range of 1 – 20 pulses per second, the **Parapulser®** provides an unrivalled range of intensities and effects.
- **0.66 Tesla pulse strength** : Each magnetic pulse has a strength of 0.66 Tesla (6,600 Gauss) at the output-surface regardless of the pulse frequency.
- **5000 magnetic pulses** : Using state-of-the-art technology, the **Parapulser®** generates up to 5,000 magnetic pulses at a time, even at the highest pulse frequency, before automatically switching off due to self-heating and thus protecting against skin burns.
- **Acoustic pulse counter** : The acoustic pulse counter enables precise and equal dosing of the magnetic pulses over the desired body locations. The acoustic pulse counter can be programmed between 10 and 9999 pulses, whereby after each counting interval a beep tone is generated.
- **Adjustable magnetic pole** : With the **Parapulser®**, the polarity of the magnetic pulses can be set to "north", "south" or "alternating pole" at the touch of a button.

- **Individual pulse frequencies** : Up to 12 individual pulse frequencies can be programmed to three digits after the dot in the range of 1 – 20 Hz.
- **Single button operation** : The **Parapulser®** can be operated with a single button in daily use.

Display functions:

- **Magnetic pole** : During operation, the polarity (north/south pole) of the generated magnetic pulse is indicated by LEDs (red/blue).
- **Intensity-step** : The currently set intensity-step (1 – 12) is displayed by color LEDs in a coded way.

Device versions:

The **Parapulser®** is available in three versions (prices see price list):

- 1.) **Parapulser® BASIC (2 Hz)** - for beginners
- 2.) **Parapulser® SEMIPRO (10 Hz)** - for common use
- 3.) **Parapulser® PRO (20 Hz)** - for professional use




The **BASIC** and **SEMIPRO** model offers a *free trial* (*) to test the enormous power of the **PRO (20 Hz)** version, if **2 Hz** or **10 Hz** are not sufficient to make a change. In addition, the **BASIC** and **SEMIPRO** model can any time be upgraded to a higher version by entering a purchased upgrade-code into the device.

(*) The free trial includes **40** applications, each with **5,000** magnetic pulses at pulse frequencies up to **20 Hz**.

Scope of delivery:

The **Parapulser®** comes in an elegant case containing all necessary parts:



Application:

Boost your well-being with help of powerful magnetic pulses and regain your energy! The **Parapulser®** supports you more effectively than any other comparable device. It's very easy to use:

- 1.) Switch on the **Parapulser®** using the "on/OFF" button,
- 2.) Start the magnetic pulses using the "on/OFF" button,
- 3.) After use, switch off the **Parapulser®** using the "on/OFF" button.



Example for application of the **Parapulser®** at the knee joint.

Legal notice:

The **Parapulser®** is not a medical device and is not intended to treat, cure or prevent any illness. It is not a substitute for medical or therapeutic treatment.

Handed out by:

Mikas Elektronik Entwicklungen

Dipl.-Ing. Alexander Mikas

Bahnhofstr. 15

82340 Feldafing

Germany

Phone : +49 8157- 92 42 92

Mail : info@mikas-elektronik.de

Web : www.mikas-elektronik.com

WEEE-Reg.-Nr. DE 13571115